

rapidshape



RAP001 £9,444.44

PRO 20

PEAK PERFORMANCE AND UNMATCHED VERSATILITY

Exceptional print performance for over 15 dental workflows and validated with over 200 dental resins - The PRO 20 delivers peak precision 3D prints at high speeds with the Rapid Shape optional patented Force Feedback technology and the newly introduced RS Turbo, printing aligner models (150µm layers) in just 7 minutes and Crowns & Bridges models (50µm layers) in 15 minutes, making it the fastest DLP 3D printer for dentistry.

The newly introduced PRO 20 RS Crystal Polish Reservoir delivers high-shine for splints, eliminating the need for post-process polishing.

Paired with the PRO wash and PRO cure, this integrated system, featuring RFID checks, ensures flawless results, seamless connectivity, and effortless delegation throughout the workflow.

Print volume: (W x D x H): 133 x 75 x 115 mm

Overall Size: 334 x 357 x 544 mm

PRO 20 ACCESSORIES

RAP006	Force Feedback/RS Turbo	£1,250.00
RAP008	PRO 20 Resin Reservoir x20	£268.52
RAP009	PRO 20 Build Platform x20	£388.89

PRO 20 BUNDLES I & II

RAP025B	Pro 20/ PROwash/ PROcure/ Accessories/ CARE Plan	£14,768.00
RAP026B	Pro 20/ PROwash/ PROcure/ Accessories/ CARE Plus	£17,546.00



RAP002 £14,074.07

D30+

SEMI-AUTOMATED PRECISION 3D PRINTER

The D30+ is a semi-automated workhorse for continuous 3D printing in dentistry, designed for unmatched efficiency and performance. Featuring the optional patented Automatic Separation Module (ASM), it automatically queues, separates, and starts the next print job, minimizing manual intervention and maximizing productivity. With integrated Force Feedback and built-in resin heating, the D30+ delivers rapid, flawless prints in any environment. Validated with over 200 resins, it offers complete material freedom with guaranteed quality. Built to excel and endure, the D30+ is your ultimate workhorse for non-stop, high-precision dental 3D printing.

Clean and cure your parts directly after printing with the PRO wash and PRO cure post-processing units for an optimal and validated end-result.

Print volume: (W x D x H): 133 x 75 x 155 (with ASM: 70 mm)

Overall Size: 480 x 410 x 690 mm

D30+ ACCESSORIES

RAP006	Force Feedback/RS Turbo	£1,250.00
RAP014	D30+ Resin Reservoir x30+	£388.89
RAP015	D30+ Build Platform x30+	£398.15

D30+ BUNDLES I & II

RAP027B	D30+/ PROwash/ PROcure/ Accessories/ CARE Plan	£18,472.00
RAP028B	D30+/ PROwash/ PROcure/ Accessories/ CARE Plus	£21,250.00